

Tests of CP and CPT symmetry with positronium

T. Namba,

K. Nishihara, T. Yamazaki, S. Asai, and T. Kobayashi
ICEPP, University of Tokyo

This talk is mainly focused on the CP test

Tests of CP and CPT symmetry
with positronium

T. Namba,

K. Nishihara, T. Yamazaki, S. Asai, and T. Kobayashi

ICEPP, University of Tokyo


Introduction

CP asymmetry

- CP asymmetry exists in the **quark sector**: $\sim 10^{-3}$ level
- How about the **lepton sector**?
- No asymmetry discovered,
but its uncertainty is $\sim 10^{-2}$ level
- Need more precision test!
- Purpose of this experiment:
Test CP symmetry with **positronium** at the 10^{-3} level

CP test with Positronium

- Angular correlations in o-Ps 3- γ decay:




CP odd term! $\left(\hat{s} \cdot \hat{k}_1 \right) \left(\hat{s} \cdot \hat{k}_1 \times \hat{k}_2 \right) \neq 0$ if CP violation

Need to measure the spin and γ -rays precisely


The simplest setup

- Spin aligned o-Ps source
- Two γ -ray detectors



- Count rate difference between N_+ and N_- means CP asymmetry

Our detector design




Magnet for the spin alignment

- Using Neodymium magnets (0.35T)


- Tensor polarization $(\hat{s} \cdot \hat{k}_1)(\hat{s} \cdot \hat{k}_1 \times \hat{k}_2)$
(spin term appears twice in the angular correlation)

- Select $m_z = \pm 1$ component
(Eliminate $m_z = 0$ component)



ortho-para mixing under the magnetic field


- $m_z=0$ component of o-Ps mixes with p-Ps
 ➡ $m_z=0$ has shorter Lifetime
- $m_z=\pm 1$ component lifetime: no change




Make delayed
timing window

↓
 $\left\{ \begin{array}{l} m_z=0: \text{eliminated} \\ m_z=\pm 1: \text{selected} \end{array} \right.$

o-Ps decay curve in a test setup (@ 0.22T)




test magnet (0.22T)





Better separation will be
obtained in the final
setup (0.35T)
 $m_z=0$ contam. < 10%

o-Ps creation & tagging system

- β^+ source: ^{22}Na ($E_{\text{end}}=546\text{keV}$), 1MBq
- β^+ stopper & e^- supplier: Silica aerogel ($\rho=0.1$)
- Tag for the β^+ emission: $t=0.1\text{mm}$ plastic scintillator



Observed spectrum with the tagging system



- Enough photoelectron collection for β^+ peak (25keV $\sim$ 11 p.e.)
- Timing resolution: $\sigma \sim 0.85\text{ns}$ (for all coincided events)

γ -ray detectors




- Requirements for the γ -ray detector:

- Large stopping power for $\sim 500\text{keV}$ γ 's
- Good energy resolution
- Good timing resolution



- LYSO ($\text{Lu}_{1.8}\text{Y}_{0.2}\text{SiO}_5$) scintillators are used for the detectors



- Large Z, high density (Atomic number of Lutetium=71, $\rho=7.1$)
- $\Delta E=10\%$ (FWHM) at 662keV
- Short decay time 40ns
- $\Delta T=150\text{ps}$ (σ) at 511keV




Turntable

- To reduce systematic errors, N_+/N_- should be compared with the same detector
(10 μm size difference can make critical asymmetry!)
- Rotating magnetic field can change the CP violation term

rotating around the axis



In reality,
we rotate the detectors, not the magnet



Current status and future prospects

- o-Ps source & β^+ -tagging system:
 - performance test was done & ready
- γ -ray detectors:
 - performance test was done & ready
- Magnet & return yoke:
 - Prototype (0.22T) was built and checked
 - Final design will be finished soon
- Turntable:
 - Final design was fixed, and now constructing the table
 - Motor control system is ready

DAQ will be started from this spring.

Estimated o-Ps event rate in the selection window $\sim 10\text{Hz}$

10^8 events will be taken in 100 days

CP test at the 10^{-3} level will be done in this year

CPT test with a similar setup

- *CPT* violating angular correlations: $\hat{s} \cdot (\hat{k}_1 \times \hat{k}_2)$
- Spin alignment (vector polarization):
 - Positrons from β^+ source are naturally polarized
 - $\langle P \rangle \sim 90\%$ (^{68}Ga source)
- For the γ -ray detection, same techniques can be used
 - LYSO scintillator pairs
 - Turntable for the systematic error reduction
- Will be measured after the *CP* test finished

Summary

- New experimental design to test CP symmetry
 - Magnet for tensor polarization
 - Plastic scintillator foils for β^+ tagging
 - LYSO crystals for γ detection
 - Turntable to reduce systematic errors
- Detector construction is almost done, and the 10^{-3} level test will be done in this year
- CPT can also be tested with similar setup